

Medical Surgical Line



Wide range of common disposables
Full CE and/or ISO certifications
Also registered in USA and in many countries
Standard languages: English-Spanish
Latest technology, materials and specifications
Produced and tested for highest quality
Solid and safe packaging with uniform appearance
Safety and confidence

Table of Contents

ADMINISTRATION OF SOLUTIONS

1. Macrodrop I.V. administration sets (20 drops/ml)
Microdrop I.V. administration sets (60 drops/ml)
Burette
Transfusion set
2. Syringe 3-60ml
Syringe 200ml in kits for contrast media
3. Syringe 0.5-1ml
Syringe 60ml, catheter tip
4. Hypodermic needles
Spinal needles
5. Stopcocks, 2 and 3-ways
I.V. catheters

WOUND CARE

6. Non-woven gauze
Gauze
7. Net bandage
Elastic bandage
8. Adhesive tapes
Nubederm adhesive Polyurethane tape

IMMOBILIZATION

9. Plaster of Paris bandage (POP)
Nucast Fiberglass bandage
Cast padding
Stockinette

STERILIZATION PAPER

10. STERIPAP flat, gusseted rolls
STERIPAP pouches

PROTECTION PRODUCTS

11. Sterile surgical and exam gloves
Non-sterile exam gloves
Caps, face masks and shoe covers

URINE DRAINAGE AND COLLECTION SYSTEM

12. Foley catheter, latex and 100% silicone
Urine bags 2000ml and 100ml

SOFT CATHETERS AND TUBES

13. Feeding tubes
Stomach Levin tubes
Urethral Nelaton catheter
Suction catheter

RESPIRATORY PRODUCTS

14. Endotracheal tube
Nebulizer mask
15. Airway devices:
Oxygen nasal canula
Oxygen mask
Oxygen mask with Venturi
Volume spacer

MORE DISPOSABLES

16. Surgical blade
I.D. bracelet
Condoms
Umbilical clamps
Digital thermometer
17. Electrodes
Dental line:
Dental needle
NuBib patient bib
Steripap sterilization pouch
Saliva aspirator
Dental gauze

Continuous and precise infusion of liquids and parenteral solutions

Macrodrop set



REFERENCE	ADDITIONAL CHARACTERISTICS
N – AAACQ-LF	Macrodrop I.V. set 20drops/ml, bulb, Y-port, luer-slip with needle 21G x 1 1/2" - latex-free
N – ABC-LF	Macrodrop I.V. set 20drops/ml, Y-port, luer-lock - latex-free
N – AWC	Macrodrop I.V. set 20drops/ml, solution filter, air inlet, Y-port, luer-slip/luer lock
N – AA2C	Macrodrop I.V. set 20drops/ml, bulb, Y-port, luer-slip
N – AVC	Macrodrop I.V. set 20drops/ml, luer-slip

1

Microdrop set



REFERENCE	ADDITIONAL CHARACTERISTICS
N – AA4C	Microdrop I.V. Set, 60 drops/ml, Y-port, bulb, luer slip
N – AALC	Microdrop I.V. Set, 60 drops/ml, solution filter, air inlet, Y-port, luer lock/luer slip

- Strong and sharp spike
- 60 drops/ml, 20 drops/ml
- Standard tube 180cm
- Solid connections
- Quantitative roller clamp for precise drip
- Y-port, bulb and needle optional
- Luer slip, luer lock or combination connector
- Latex-free available
- Air inlet and/or solution filters available
- Cellophane packaging, easy to open
- 5 year expiry date

CE 0197

Maximum dosage precision and volumetric administration of liquids and medication

Burette



- Hypo-allergenic silicone floater available
- Soft chamber with white strip for writing and visualization
- Blue graduated volumetric scale 0 to 150ml
- Hydrophobic and bacteriostatic filter
- Strong Nitrile spike
- Hanger for greater stability
- Latex-free injection port
- Label with complete instructions
- Cellophane packaging, easy to open
- 5 year expiry date

CE 0197

REFERENCE	ADDITIONAL CHARACTERISTICS
N – AA9F - LF	Burette, 20 drops/ml, for connection to I.V. set
N – AQF - LF	Burette, 20 drops/ml, for connection to I.V. set with floater
N – AADD - LF	Burette, complete set, 60 drops/ml, LL/LS combined connector
N – AHD-LF	Burette, complete set, 60 drops/ml, luer slip connector

Transfusion set



- Designed to customer specification
- Strong and sharp spike
- 15 and 20 drops/ml available
- Double drip soft chamber
- Top chamber with a 200 µg microfilter
- Standard 180cm tube, little memory
- Quantitative roller clamp for precise drip
- Latex-free
- Y-port, air inlet, bulb and needle optional
- Luer-slip, luer-lock or combination connector
- Cellophane packaging, easy to open
- 5 year expiry date

REFERENCE	ADDITIONAL CHARACTERISTICS
N – AACB - LF	20 drop/ml transfusion set. Luer slip, latex-free

Nubenco®

Administration of solutions

SYRINGES 3 - 60ml



Ultra accuracy in dosage and safety

2

- Our syringes have 3 pcs: barrel, plunger, rubber stopper
- Attached needle and latex-free optional
- Luer-lock or luer-slip
- Breakable plunger - NO RE-USE
- Siliconized for easy and smooth plunger slide
- Positive plunger stop
- Permanent engraved fine scale marking on barrel
- Grid plunger top for controlled injection
- 5 year expiry date

CE 0197

REFERENCE	DESCRIPTION
N - S4H	3 ml, luer slip, no needle
N - S1HQ	3 ml, luer lock, needle 21G x 1 1/2"
N - S1HT	3 ml, luer lock, needle 22G x 1 1/2"
N - S1HV	3 ml, luer lock, needle 23G x 1"
N - S4I	5 ml, luer slip, no needle
N - S1IW	5 ml, luer lock, needle 23G x 1 1/2"
N - S1IQ	5 ml, luer lock, needle 21G x 1 1/2"
N - S4J	10 ml, luer slip, no needle
N - S1JQ	10 ml, luer lock, needle 21G x 1 1/2"
N - S1KQ	20 ml, luer lock, needle 21G x 1 1/2"
N - S3L	20 ml, luer lock, no needle
N - S3M	30 ml, luer lock, no needle
N - S3P	50 ml, luer lock, no needle
N - S4P	50 ml, luer slip, no needle
N - S3R	60 ml, luer lock, no needle
N - S4R	60 ml, luer slip, no needle

Nubenco is a registered trademark of Nubenco Enterprises Inc. - Address: 1 Kalisa Way, Paramus, New Jersey 07652 - USA. Ph.: (201) 967-9000 - www.nubenco.com - E-mail: info@nubenco.com

Nubenco®

Administration of solutions

SYRINGE 200ml FOR CONTRAST MEDIA

- High quality Polypropylene material
- For administration of Contrast media
- Kit for CT injector systems:
 - Syringe 200ml
 - Spiral connector tube
 - Filling tube
 - Detailed insert
 - Individually boxed
- Fits perfectly to CT injectors such as Medrad 200/300/Vistron series and Stellant models*
- Connector tube with female/male adapter; 150cm length available separately (Ref. N-Q2)
- Double packaging for better protection and storage
- Stands pressures up to 355 PSI/2446 kPa
- 5 year expiry date



CE 0197



Nubenco is a registered trademark of Nubenco Enterprises Inc. * Stellant, Vistron, Medrad are registered trademarks of Medrad Address: 1 Kalisa Way, Paramus, New Jersey 07652 - USA. Ph.: (201) 967-9000 - www.nubenco.com - E-mail: info@nubenco.com

Nubenco®

Administration of solutions

Nubenco®

Specialized needles and syringes

Syringe 0.5 - 1ml

Ultra accuracy in dosage and safety



- Our syringes have 3 pcs: barrel, plunger and rubber stopper
- Luer slip connector and detachable needle
- For 0.5ml and 1ml, fixed needle available for high precision of insulin and tuberculin dosage
- All gauges, lengths - color-coded
- Sharp bevel and siliconized needle for less traumatic puncture
- Longer tip for minimal loss during injection
- Siliconized for easy, smooth plunger slide
- Positive plunger stop
- Permanent engraved markings on barrel
- Plunger flat top with grids for controlled injection

3

REFERENCE	GAUGE	LENGTH	COLOR	APPLICATION
N-S2F5	30G	1/2"	Yellow	For Insulin, cosmetic use (Botox)
N-S2FAI	29G	1/2"	Red	
N-S2F4	28G	1/2"	Blue	
N-S2F3	27G	1/2"	Light Gray	
N-S2F2	26G	1/2"	Brown	Allergic test, mesotherapy Intradermic use, subcutaneous
N-S2F1	25G	5/8"	Orange	
N-S9ZAZ	30G	5/16"	Beige	Vaccine Fixed needle Intradermic and subcutaneous use
N-S2TV	23G	1"	Dark Blue	
N-S2TAC-LF	22G	1 1/4"	Gray	

CE 0197

For feeding, irrigation and internal wash


- 60ml syringe
- Universal tip to fit different gauges tubes or catheters
- Includes separate luer slip adapter
- Siliconized for easy and smooth plunger slide
- Positive plunger stop
- Grid plunger top for controlled irrigation

Catheter tip syringe 60ml



REFERENCE	ADDITIONAL CHARACTERISTICS
N - SDR - LF	Syringe 60ml, catheter tip, luer slip adapter

Nubenco is a registered trademark of Nubenco Enterprises Inc. - Address: 1 Kalisa Way, Paramus, New Jersey 07652 - USA. Ph.: (201) 967-9000 - www.nubenco.com - E-mail: info@nubenco.com



**Painless
Better quality of life
and economy**

Nubenco Fixed needle "Zero dead space" (0.015cc Loss)
Nubenco Detachable needle with long rubber tip (0.04cc Loss)
Normal Detachable needle with regular rubber tip (0.06cc Loss)

* Wasted solution left in syringe after administration

- Available for insulin, immunization, cosmetology, allergy, tuberculin, etc
- "Painless" needles, sharp tip and siliconized
- Siliconized rubber tip and barrel slide smoothly
- Long rubber tip to reduce wasted solutions with detachable needle
- Maximum cost savings with fixed needle for insulin and tuberculin, 0.5 - 1ml
- Nubenco 29G-30G needles fill 1ml in just 15 seconds
- Insulin markings of 100IU, 80IU, 50IU or 40IU
- Insulin syringes with detachable needle 25G-30G
- Also insulin and tuberculin syringes with fixed needle 27G-30G
- 5 year expiry date

Nubenco syringes with detachable needle have a rubber long "special extended" tip that fits perfectly into the needle hub fitting to reduce the liquid lost in this fitting. Nubenco syringe with fixed needle has a "zero dead space" with minimum losses of liquid in the hub of needle. For both needles, the filled volume is measured by the volumetric markings at the superior section of the rubber tip, without checking the extension of tip

Nubenco is a registered trademark of Nubenco Enterprises Inc. - Address: 1 Kalisa Way, Paramus, New Jersey 07652 - USA. Ph.: (201) 967-9000 - www.nubenco.com - E-mail: info@nubenco.com

HYPODERMIC needle
Easy puncture, minimal pain

Ultra accuracy in dosage and safety



4

CE 0197

REFERENCE	COLOR	DESCRIPTION
N - N1B	White	Needle 16G x 1/2"
N - N1A	White	Needle 16G x 1"
N - N19	White	Needle 16G x 1 1/2"
N - N1F	Pink	Needle 18G x 1"
N - N1G	Pink	Needle 18G x 1 1/2"
N - N1I	Beige	Needle 19G x 1"
N - N1J	Beige	Needle 19G x 1 1/2"
N - N1L	Yellow	Needle 20G x 1"
N - N1M	Yellow	Needle 20G x 1 1/2"
N - N1P	Green	Needle 21G x 1"
N - N1Q	Green	Needle 21G x 1 1/2"
N - N1S	Black	Needle 22G x 1"
N - N1T	Black	Needle 22G x 1 1/2"
N - N1V	Blue	Needle 23G x 1"
N - N1W	Blue	Needle 23G x 1 1/2"
N - N1Y	Purple	Needle 24G x 1"
N - N1AE	Orange	Needle 25G x 1"
N - N1I	Orange	Needle 25G x 3/8"
N - N12	Brown	Needle 26G x 1/2"
N - N13	Gray	Needle 27G x 1/2"
N - N1AI	Salmon	Needle 29G x 1/2"
N - N15	Ice	Needle 30G x 1/2"

- Fits luer slip and luer lock connector
- High grade stainless steel, "painless"
- Sharp bevel and siliconized cannula
- Various gauges 16-30G
- Color coded hub
- Barcode printed on blister
- Sterile individual blister
- 5 year expiry date

SPINAL NEEDLES

- For the administration of anesthesia in the spinal area
- High grade stainless steel cannula
- Different types:
 - Whitacre - pen tip, atraumatic, with stylet introducer, for intradural blockage and LCR sampling
 - Tuohy - Plastic cannula with insertion markings (each 1cm) to control insertion depth
- Sterile blister packaging "easy open", rigid protector cap
- 5 year expiry date

Premium Spinal and Epidural needles



CE 0197

2 and 3-way Stopcocks



5

CE 0197

- Available in white, red or blue color
- Rotational knob with flow arrows (360°)
- Connects to all luer-slip and luer-lock devices
- With extension sets, any length available
- Produced with medical grade PVC, transparent
- High resistance to mechanical or manual flow pressure
- Sterile individual blister packaging
- 5 year expiry date


REFERENCE	DESCRIPTION
N - QTAE	3-way stopcock with extension, 2.6mm x 15cm
N - QTAF	3-way stopcock with extension, 2.6mm x 50cm
N - QTAG	3-way stopcock with extension, 2.6mm x 100cm
N - QUAJA	3-way stopcock, 2.6mm, color red
N - QUAJB	3-way stopcock, 2.6mm, color blue

- Sterile, single-use
- Available in pen-type or blister
- Available wings, air inlet and injection port
- Available with heparin cap and connector for closed infusion
- Available with large or small needle cap
- Available new Y-type
- High grade stainless steel needle
- X-ray detectable plastic needle
- Needle gauges from 14G to 24G
- 5 year expiry date

EASY INSERTION PAINLESS



Y-type




Injection valve type



Pen-type model

CE 0197

- Economic and very convenient for routine applications
- Ideal for wound care: Protection, cleansing and dressing
- Skin-friendly; low adhesion to tissues
- Highly absorbent
- Hypoallergenic
- Available also with drain for Tracheostomy and Phlebotomy
- Sterile 4-10pcs per blister or 2-3pcs per envelope
- Sterile - 5 year expiry date




REGULAR: 5 x 5cm (2" x 2")
7.5 x 7.5cm (3" x 3")
10 x 10cm (4" x 4")

Rayon

Polyester

30 or 40 g/m²



WITH DRAIN: 5 x 5cm (2" x 2")
10 x 10cm (4" x 4")

CE 0197

- Excellent quality for wound care
- Good air circulation and high absorbent capacity
- Folded and unfolded edges
- USP Type VII, varying weights for routine use
- USP Type V for greater absorption
- Available 4, 6, 8, 12 and 16ply
- 4"x4" X-ray detectable, for open surgery - in blister packaging of 5-10pcs
- Rolls 36" x 90m : USP type VII - 2 ply
USP type V - 4 ply
- Rolls 1", 2", 3", 4", 6" x 5 yds, non-sterile, medical grade paper packaging
- Non-woven gauze, 30-40 g/m 4 and 6ply
- Sterile - 5 year expiry date

All types of natural gauzes,
sterile and non-sterile



CE 0197

Replaces adhesive tapes for greater patient freedom, mobility and user convenience

- Breathable cover for gauze dressings
- Washable and reusable
- Width to fit every body part
- Can be cut and shaped
- Lightweight
- Stretchable
- Net bandage cap available, Ref. N-G653



7

100% Cotton
Tubular and hypoallergenic



All widths available:

- 0.8cm x 10yds, #1
- 1.6cm x 10yds, #2
- 2.7cm x 10yds, #3
- 3.5cm x 10yds, #4
- 5.1cm x 10yds, #5
- 6.2cm x 10yds, #6
- 7.5cm x 10yds, #7

**Complete range of elastic bandages
Wide variety of sizes
High quality**

- Elastic Bandage Plain
- Elastic Bandage Crepe
- Supplied with two metal clamps



- ½", 1", 2", 3", 4", 6" x 1, 5, 10yds
- Easy to tear or cut
- Hypoallergenic
- Excellent adhesion
- Skin friendly
- Porous and breathable
- Dispensers, spools or tubes



8

Plastic - NUPLAST®
Paper - NUPAP®
Silk - NUSEDA®
Zinc Oxide - NUTELA®
Non-woven dressing - NUFIX®



CE 0197

POLYURETHANE ADHESIVE TAPE

- Water permeable, breathable, conformable
- Accelerates healing, helps to prevent infection
- Long lasting adhesion
- Hypoallergenic adhesive
- Available shapes: Rectangle, square and oval
- Supplied with or without central dressing
- Size 6x7cm available with IV catheter slot
- 3 year expiry date

Feels likes a second skin



CE 0197

Full line of quality and reliability

- Complete immobilization line
- Size: 2", 3", 4", 5", 6", 8" x 3-4-5 yards
- POP setting time: 5-8 minutes
- Fiberglass setting time: 3-5 minutes
- Synthetic Stockinette (Terylene / Cotton), for lower absorption, greater dryness and better air permeability, preferred over 100% cotton
- Acrylic Cast Padding (Rayon), acts as an air permeable cushion system controlling temperature changes, preferred over other materials

NUCAST®
Fiberglass bandage



Stockinette



Cast padding



POP bandage

9

P.O.P. PLASTER OF PARIS BANDAGE

- Fast setting time 5-8 Minutes
- Drying time 24 Hours, for greater patient comfort
- Variety of sizes
- Provides a strong and safe immobilization
- High resistance, avoiding damage and increasing patient comfort
- Lower plaster losses, for higher resistance and better consistency
- High creaminess, providing a superior finish
- Radiolucent for X-ray imaging and diagnostics
- Air-tight packaging, preserving the product from humidity

Excellent quality and reliability

Dimensions:

3" x 5 yds


4" x 5 yds

5" x 5 yds

6" x 5 yds

8" x 5 yds





10

**Indicators for vapor and ethylene oxide gas. One side medical grade paper and other side transparent film polyester-polyethylene.
Flat or gusseted rolls and pouches in all available sizes**

- Easy to use, clear view of products inside
- Longer sterility period
- 60g/m² high density paper
- Indicators along edges of paper roll
- Repels liquids
- Antimicrobial barrier
- Cuts and seals reliably
- No ETO residue
- Complies with EN 868

Medical Grade Paper



Available sizes and types

REFERENCE	SIZE (metric system)	ROLL TYPE	PACKAGING
N-PF1A	50mm x 200m	FLAT	6 rolls/box
N-PF1B	75mm x 200m	FLAT	4 rolls/box
N-PF1C	100mm x 200m	FLAT	3 rolls/box
N-PF1D	150mm x 200m	FLAT	2 rolls/box
N-PF1E	200mm x 200m	FLAT	1 rolls/box
N-PF1F	250mm x 200m	FLAT	1 rolls/box
N-PF1G	300mm x 200m	FLAT	1 rolls/box
N-PF1H	350mm x 200m	FLAT	1 rolls/box
N-PF1I	400mm x 200m	FLAT	1 rolls/box

REFERENCE	SIZE (metric system)	ROLL TYPE	PACKAGING
N-PG1C	100mm x 100m	GUSSETED	2 rolls/box
N-PG1D	150mm x 100m	GUSSETED	1 rolls/box
N-PG1E	200mm x 100m	GUSSETED	1 rolls/box
N-PG1F	250mm x 100m	GUSSETED	1 rolls/box
N-PG1G	300mm x 100m	GUSSETED	1 rolls/box
N-PG1H	350mm x 100m	GUSSETED	1 rolls/box
N-PG1I	400mm x 100m	GUSSETED	1 rolls/box

LATEX GLOVES

Sterile gloves



For protection of medical staff and patients in aseptic procedures

- 100% latex
- Low protein content
- High resistance and elasticity
- Textured or plain
- High sensitivity
- Lightly powdered with corn starch
- Powder-free available

11

REFERENCE	PRODUCT	DESCRIPTION
N-L1A1	Sterile glove 6.0	Box x 50pairs
N-L1B1	Sterile glove 6.5	Box x 50pairs
N-L1C1	Sterile glove 7.0	Box x 50pairs
N-L1D1	Sterile glove 7.5	Box x 50pairs
N-L1E1	Sterile glove 8.0	Box x 50pairs
N-L1F1	Sterile glove 8.5	Box x 50pairs
N-L2A	Sterile exam glove 6.0	Box x 50pairs
N-L2B	Sterile exam glove 6.5	Box x 50pairs
N-L2C	Sterile exam glove 7.0	Box x 50pairs
N-L2D	Sterile exam glove 7.5	Box x 50pairs
N-L2E	Sterile exam glove 8.0	Box x 50pairs
N-L2F	Sterile exam glove 8.5	Box x 50pairs

CE 0197

Exam Gloves

Extra protection for medical staff in routine medical procedures, and for industrial applications

STERILE:

- Surgeon's: left and right
- Exam: ambidextrous
- Double packaging per pair
- 3 year expiry date

NON-STERILE:

- Ambidextrous
- With independent certification
- PE, Vinyl and Nitrile available



REFERENCE	PRODUCT	DESCRIPTION
N - L3G	EXAM GLOVE, NON-STERILE, SIZE S	BOX OF 100
N - L3H	EXAM GLOVE, NON-STERILE, SIZE M	BOX OF 100
N - L3I	EXAM GLOVE, NON-STERILE, SIZE L	BOX OF 100

Caps, Face mask and Shoe cover

Protection gear that blocks particulate and biological contamination in Medical and Industrial areas



- Wide range of applications
- Hypoallergenic and hydrorepellant polypropylene material, non woven of excellent quality
- Smooth and anatomical fit
- Available in colors: green, blue, white
- Even texture, dissipates heat
- Disposable surgical caps with ties
- Elastic cap for Nurses and other specialized staff
- Face masks with 3 layer biological filter, with ties or earloops and bendable nose strip
- Slip-free shoe cover, adjusts to footwear

REFERENCE	PRODUCT	DESCRIPTION
N - CD2	Surgeon cap	BOX x 100
N - CD1	Surgeon cap	BOX x 100
N - CG2	Cap for nurse or surgical technician, round	BOX x 100
N - CV1	Accordeon nurse cap	BOX x 100
N - CA4	Face mask, with ties	BOX x 100
N - CB4	Face mask, with ear loops	BOX x 100
N - CG4	Shoe cover	BOX x 100

Foley catheter with balloon

CE 0197



Great comfort for urinary drainage and irrigation

Wide range of catheters for therapeutic drainage, suction, infusion

- Latex coated in silicone, for smooth insertion and comfort for user
- Optimum balance of hardness and softness
- Reinforced juncture for greater security
- 2-way in gauges 6, 8, 10, 12, 14, 16, 18, 20, 22 and 24FR
- 3-way in gauges 18, 20, 22, 24, 26FR
- 6, 8, 10FR with rigid catheter guide
- Our 3-way catheters are excellent for bladder irrigation
- Fits luer slip and luer lock connector, all collection and catheter tip syringe systems
- Standard length 40cm, 2 recessed side openings at balloon tip
- Color coded inflating valve
- Symmetric balloon 3-50ml volumes available
- Blister sterile packaging with inner liner - 5 year expiry date
- 100% Silicone Foley catheter available, all sizes and gauges

12

LATEX REFERENCE	SILICONE REFERENCE	COLOR	GAUGE	DESCRIPTION
N-QQH1		BLACK	6 FR	2 way, balloon 3-5ml, 30cm
N-QQL1		BLACK GRAY	8 FR	2 way, balloon 3-5ml, 30cm
N-QQQ1		LIGHT GRAY	10 FR	2 way, balloon 3-5ml, 30cm
N-QQR2		WHITE	12 FR	2 way, balloon 5-15ml, 40cm
N-QOS2		LIGHT GREEN	14 FR	2 way, balloon 5-15ml, 40cm
N-QOT2	N-Q7T3	ORANGE	16 FR	2 way, balloon 5-15ml, 40cm
N-QOU2	N-Q7U3	RED	18 FR	2 way, balloon 5-15ml, 40cm
N-QOV2		YELLOW	20 FR	2 way, balloon 5-15ml, 40cm
N-QOW2		PURPLE	22 FR	2 way, balloon 5-15ml, 40cm
N-QOX2		BLUE	24 FR	2 way, balloon 5-15ml, 40cm
N-QWU4		RED	18 FR	3 way, balloon 30ml, 40cm
N-QWV4		YELLOW	20 FR	3 way, balloon 30ml, 40cm
N-QWW4	N-Q8W5	PURPLE	22 FR	3 way, balloon 30ml, 40cm
N-QWX4	N-Q8X5	BLUE	24 FR	3 way, balloon 30ml, 40cm
N-QWY4		BLUE	26 FR	3 way, balloon 30ml, 40cm

CE 0197

Specialized product



N-U9

N-U9

- Most sophisticated and thicker bag
- Ideal for urologic post-surgery
- 1cm diameter tube to reduce risk of clogging
- 110cm long tube
- Convenient sampling port
- Double hanging support-chord and velcro
- Side pockets to protect drainage tube
- N-U9-NR : Bag without anti-reflux valve, faster flow, great for ER applications

N-U7 & N-U10

- 0.5cm diameter tube
- 95cm long tube
- N-U10, clip-on hanging device

100ml bag size for safe pediatric procedures



N-U3

Comfort, safety and precision for drainage and urinary quantification

- Different bags for various conditions
- Strong 2000 ml bag, medical grade PVC
- Supplied with safe hanging system
- Universal connector fits all catheters
- Soft entry tube for greater patient comfort
- NRV - non-reflux valve available
- Transparent bag front for easy volume inspection
- Precise volumetric markings, bright blue
- Patient information area
- Easy, reliable fixed drainage valve with protector cap for greater hygiene
- Sterile available - 5 year expiry date
- 500ml leg bag with strap available

Safety and economy in urinary drainage systems



N-U7



N-U10

N-U3

- Sterile, transparent PE material
- 10-100ml volume capacity
- Visible blue markings
- Hypoallergenic adhesive and protective pad

Feeding tube

Wide range of applications such as infusion, drainage and other therapeutic objectives

For enteral nutrition and other gastric procedures

- Luer slip connector with protector cap
- 6, 8 and 10FR gauges available
- Soft PVC with open tip and 2 side openings
- 5 year expiry date

REFERENCE	COLOR	DESCRIPTION
N-QFH	DARK BLUE	Feeding tube, 6FR, 40cm
N-QFL	LIGHT BLUE	Feeding tube, 8FR, 40cm
N-QFQ	CLEAR	Feeding tube, 10FR, 120cm

CE 0197

13

Urethral Nelaton catheter

For drainage, urine sampling, bladder decompression

- Medical grade PVC, soft and transparent
- Universal conical connector
- 5, 6, 8, 10, 12, 14, 16, 18, 20, 22FR x 40cm
- 2 side openings at tip
- 5 year expiry date

REFERENCE	COLOR	DESCRIPTION
N - QKF	YELLOW	Nelaton catheter, 5 FR, 40cm
N - QKH	LIGHT GREEN	Nelaton catheter, 6 FR, 40cm
N - QKL	BLUE	Nelaton catheter, 8 FR, 40cm
N - QKQ	BLACK	Nelaton catheter, 10 FR, 40cm
N - QKR	WHITE	Nelaton catheter, 12 FR, 40cm
N - QKS	DARK GREEN	Nelaton catheter, 14 FR, 40cm
N - QKT	ORANGE	Nelaton catheter, 16 FR, 40cm
N - QKU	RED	Nelaton catheter, 18 FR, 40cm
N - QKV	YELLOW	Nelaton catheter, 20 FR, 40cm
N - QKW	DARK BLUE	Nelaton catheter, 22 FR, 40cm

CE 0197

Allows gastric decompression, sampling, wash, drainage and feeding

Stomach Levin tube

- Universal conical adapter
- Flexible tube with open tip and 4 side openings
- Gauges: 6,8,10,12FR x 90cm
- 10,14, 16, 18, 20FR x 90cm or 120cm
- Markings on tube for accurate insertion
- X-ray line optional (Cat. REF N-QI*)
- 5 year expiry date

REFERENCE	COLOR	DESCRIPTION
N - QACH	LIGHT GREEN	Stomach Levin tube 6FR x 90cm
N - QHL	BLUE	Stomach Levin tube 8FR x 125cm
N - QHQ	BLACK	Stomach Levin tube 10FR x 125cm
N - QHR	WHITE	Stomach Levin tube 12FR x 125cm
N - QHS	GREEN	Stomach Levin tube 14FR x 125cm
N - QHT	ORANGE	Stomach Levin tube 16FR x 125cm
N - QHU	RED	Stomach Levin tube 18FR x 125cm
N - QHV	YELLOW	Stomach Levin tube 20FR x 125cm

Eases controlled suction of bodily fluids and secretions

Suction catheter

- Thumb-controlled connector, fits any suction system
- Easy gauge identification - color coded
- Gauges: 5, 6, 8, 10, 12, 14, 16, 18FR x 60cm
- Medical grade PVC, transparent
- Open tip with 2 side openings
- Markings for controlled insertion
- 5 year expiry date

REFERENCE	COLOR	DESCRIPTION
N - QEF	YELLOW	Suction catheter 5FR, 50cm
N - QEH	GREEN	Suction catheter 6FR, 50cm
N - QEL	BLUE	Suction catheter 8FR, 50cm
N - QEQ	BLACK	Suction catheter 10FR, 50cm
N - QER	WHITE	Suction catheter 12FR, 50cm
N - QES	DARK GREEN	Suction catheter 14FR, 50cm
N - QET	ORANGE	Suction catheter 16FR, 50cm
N - QUE	RED	Suction catheter 18FR, 50cm

Assures good permeability and ventilation of the airway passage

Endotracheal tube

- PVC tube with visible markings and radiopaque line
- Different gauges with or without balloon (2-30ml capacity)
- Flexible material that adapts well in trachea
- Softens with body heat
- Murphy type
- Perfect fit to anesthesia and ventilation circuits
- 5 year expiry date

Endotracheal tube reinforced

- Stainless steel spiral wire

14

CE 0197

Device with optimum mist action for respiratory procedures

- Pediatric and adult sizes
- Fitted mask, made with soft PVC
- Large air orifices at sides of mask for CO₂ release
- Adjustable elastic band, latex-free
- Nebulizer cup, 7ml capacity, vaporizes in all positions
- 2m tube for better patient mobility

Nebulizer mask



REFERENCE	ADDITIONAL CHARACTERISTICS
N - Q49	Nebulizer mask, pediatric, Small, 200mm tube
N - Q46	Nebulizer mask, pediatric, Medium, 200mm tube
N - Q47	Nebulizer mask, adult, Large, 200mm tube
N - Q4AK	Nebulizer mask, adult, X-large, 200mm tube

Endotracheal tube with balloon

CE 0197

REFERENCE	DESCRIPTION	GAUGE
N - QCC	Endotracheal tube with balloon, PVC, radiopaque line	14 FR, 3.5mm ID
N - QCD	Endotracheal tube with balloon, PVC, radiopaque line	16 FR, 4.0mm ID
N - QCE	Endotracheal tube with balloon, PVC, radiopaque line	18 FR, 4.5mm ID
N - QCF	Endotracheal tube with balloon, PVC, radiopaque line	20 FR, 5.0mm ID
N - QCG	Endotracheal tube with balloon, PVC, radiopaque line	22 FR, 5.5mm ID
N - QCH	Endotracheal tube with balloon, PVC, radiopaque line	24 FR, 6.0mm ID
N - QCI	Endotracheal tube with balloon, PVC, radiopaque line	26 FR, 6.5mm ID
N - QCJ	Endotracheal tube with balloon, PVC, radiopaque line	28 FR, 7.0mm ID
N - QCK	Endotracheal tube with balloon, PVC, radiopaque line	30 FR, 7.5mm ID
N - QCL	Endotracheal tube with balloon, PVC, radiopaque line	32 FR, 8.0mm ID
N - QCM	Endotracheal tube with balloon, PVC, radiopaque line	34 FR, 8.5mm ID
N - QCN	Endotracheal tube with balloon, PVC, radiopaque line	36 FR, 9.0mm ID
N - QCP	Endotracheal tube with balloon, PVC, radiopaque line	38 FR, 9.5mm ID

Endotracheal tube without balloon

REFERENCE	DESCRIPTION	GAUGE
N - QDA	Endotracheal tube without balloon, PVC, radiopaque line	10 FR, 2.5mm ID
N - QDB	Endotracheal tube without balloon, PVC, radiopaque line	12 FR, 3.0mm ID
N - QDC	Endotracheal tube without balloon, PVC, radiopaque line	14 FR, 3.0mm ID
N - QDD	Endotracheal tube without balloon, PVC, radiopaque line	16 FR, 4.0mm ID
N - QDE	Endotracheal tube without balloon, PVC, radiopaque line	18 FR, 4.5mm ID
N - QDF	Endotracheal tube without balloon, PVC, radiopaque line	20 FR, 5.0mm ID
N - QDG	Endotracheal tube without balloon, PVC, radiopaque line	22 FR, 5.5mm ID
N - QDJ	Endotracheal tube without balloon, PVC, radiopaque line	28 FR, 7.0mm ID
N - QDK	Endotracheal tube without balloon, PVC, radiopaque line	30 FR, 7.5mm ID

Endotracheal tube reinforced

REFERENCE	DESCRIPTION	GAUGE
N - Q5Z	Endotracheal tube reinforced with balloon, PVC, radiopaque line	28FR, 2.5mm ID
N - Q51	Endotracheal tube reinforced with balloon, PVC, radiopaque line	30FR, 3.0mm ID
N - Q52	Endotracheal tube reinforced with balloon, PVC, radiopaque line	32FR, 3.5mm ID

OXYGEN NASAL CANNULA

Flexible and soft nosetips
Strap ring adjusts to patient permitting movement

- Four models for different situations
- Neonatal, pediatric and adult sizes
- Star lumen tube 2 m with universal connector
- Designed for maximum comfort and utility
- Soft and contoured masks, with comfortable strap and adjustable nose clip
- Latex free
- Medical grade PVC
- Individual sterile packaging

OXYGEN MASK


With open holes for exit of exhaled gas

VENTURI OXYGEN MASK

With 2 or 6 color-coded diluters to vary the O₂ concentration
With open holes for exit of exhaled gas

NEBULIZER MASK

For treatment of respiratory or lung disease
Nebulizes in all directions
Cup maximum capacity 7ml
With small holes to inhale air during treatment



CE 0197

VOLUME SPACER

- Latex-free mask comfortably fits face contour
- Adult and pediatric size, 180ml
- Prevents residual layer of inhaled medicine in lower breathing chambers
- Maintains a natural breathing rhythm
- Maximized effectiveness of medicine dosage
- One-way valve to prevent any medicine leak
- Manufactured according to essential requirements of norms EN1441: 1994 and norm ISO 11607
- Printed instructions on packaging
- Made with resistant materials, durable and reusable when thoroughly washed after use

For controlled inhalation therapy

Capacity of 180ml
One-way valve



Surgical blades



Sharp tool for precise incisions

- Carbon or stainless steel
- All sizes available
- Faster cuts with our sharp blades
- Scalpel with blade available
- Inner packaging, for protection
- Outer packaging anti-corrosive
- Sterile 3 year expiry date

- Pediatric colors: Pink and blue
- White color for adult use
- Adjustable clip

Condoms



- Lubricated latex in aluminum foil
- Reservoir tip
- Available in 3 pcs retail box, or 144pcs gross box
- Many styles/types available
- Gamma ray sterilization
- 3 year expiry date

- Designed to firmly compress the newborn's umbilical chord
- PVC, apyrogenic
- 5 year expiry date

Umbilical clamp



CE 0197

I.D. Bracelet



16


DIGITAL THERMOMETER

- Water resistant to permit easy cleaning
- Flexible tip for greater comfort
- Centigrade reading
- Faster response time: 1 minute
- Travel size and easy to use
- Fever alarm
- Last reading memory
- Auto Shutoff
- Also available: Other retail packaging
- Talking thermometer (different languages)
- With plastic sheaths



Electrodes

Excellent conductivity for perfect and precise electrocardiogram trace



17

- Hypoallergenic adhesives
- Use for patient monitoring
- Up to 72 hrs. monitoring, long durability
- Transparent gel
- No extra gel required
- Safe stainless steel clip
- Adult electrode, non-woven
- Pediatric and neonatal electrode, foam
- Excellent packaging with aluminum inside liner

REFERENCE	SIZE	PACKAGING
N-XV	Adult	100 pcs/bag
N-XY	Pediatric	100 pcs/bag
N-XZ	Neonatal	100 pcs/bag

- **Dental Gauze:**
2" x 2" 8 ply
200 pcs. Package
- **Dental Needle:**
27 G Long
30 G Short
- **Saliva Aspirator:**
Copper Wire
Rounded Multi-channeled tip
Wide diameter for better suction
- **NuBib:**
Dual Layer
Water proof
Absorbent
46cm x 58 cm
- **STERIPAP:**
Sterilization Pouches
Flat and gusseted
- **Face Mask:**
Aluminum nose clip
Elastic ear loops
Non Woven Material
3 layer filtration
Skin friendly
- **Gloves:**
Latex, Vinyl and PE
Powdered
Powdered free

High quality products for Dentistry



CE 0197