


Nubenco®

Administration of solutions and
blood components



IV CATHETER

Easy Insertion and Painless

- Sterile, single-use
- Available in pen-type or blister packaging
- Available with wings, air inlet and injection port
- Available with heparin cap and connector for closed infusion
- Available new Y-type
- High grade stainless steel needle
- X-ray detectable plastic cannula
- Needle gauges from 14G to 24G
- 5 year expiry date



Gauge	Pen-type model	Y-type	Radiopaque line, Hydrophobic filter	Injection port and wings
14G	N-I2AJ (2.0 OD, 55mm)			N-I8AJ (1.6 OD, 55 mm)
16G	N-I2AK (1.8 OD, 55mm)			N-I8AK (1.6 OD, 55mm)
18G	N-I2AL (1.3 OD, 45mm)	N-I7AL (1.3 OD, 45 mm)	N-I5B3-R (1.3 OD, 32 mm)	N-I8AL (1.2 OD, 45mm)
20G	N-I2AM (1.1 OD, 32mm)	N-I7AM (1.1 OD, 32 mm)	N-I5AL-R (1.3 OD, 45 mm)	N-I8AM (1.0 OD, 33mm)
22G	N-I2AC (0.9 OD, 32mm)	N-I7AN (0.9 OD, 25 mm)	N-I5AM-R (1.1 OD, 32 mm)	N-I8AN (0.8 OD, 25mm)
24G	N-I1D7 (0.7 OD, 33mm)	N-I7AO (0.7 OD, 19 mm)	N-I5AN-R (0.9 OD, 25 mm)	N-I8AO (0.7 OD, 19 mm)
24G	N-I1AO (0.7 OD, 19mm)		N-I5AO-R (0.7 OD, 19 mm)	

Nubenco Enterprises, Inc.

Need to contact us? Call 1.201.967.9000 or

Visit our website www.nubenco.com • e-mail: info@nubenco.com