

### Syringe 200 ml For Contrast Media



- High quality Polypropylene material
- For administration of Contrast Media.
- Kit for CT injector systems:
  - Syringe 200 ml
  - Low pressure connector tube
  - Filling tube
  - Detailed insert
  - Individually boxed
- Fits perfectly to CT injectors such as Medrad 200/300/Vistron Series and Stellant Models
- Kits also available for other systems
- Connector tube with female/male adapter; 150 cm Length available separately (Ref. N-Q2)
- Double packaging for better protection and storage.
- Stands pressures up to 355 PSI/2446 kPa
- 5 year expiry date

REFERENCE	DESCRIPTION
N - SCY	200 ml Syringe, quick fill tube, low pressure tube ( fits Medrad 200/300 series)
N - SEY	200 ml Syringe, quick fill tube, low pressure tube ( fits Stellant Type Injector )
N - SFY	200 ml Syringe, Frontal injector, quick fill tube, low pressure tube ( fits CT9000 Injector)
N - SGY	2-200 ml Syringe, quick fill tube, low pressure tube ( fits Stellant Type Injector)