

ENDOTRACHEAL TUBE



Assures good permeability and ventilation of the airway passage

- PVC tube with visible markings and radiopaque line
- All gauges available with or without balloon (2-30 ml capacity)
- Tube material softens with body heat
- Murphy type
- Perfect fit to anesthesia and ventilation circuits
- 5 year expiry date

Endotracheal tube reinforced

- Stainless steel spiral wire

GAUGE	ENDOTRACHEAL TUBE WITHOUT BALLOON	ENDOTRACHEAL TUBE WITH BALLOON	ENDOTRACHEAL TUBE REINFORCED
2.5 mm ID 10 FR	N - QDA		
3.0 mm ID 12 FR	N - QDB		
3.5 mm ID 14 FR	N - QDC	N - QCC	
4.0 mm ID 16 FR	N - QDD	N - QCD	
4.5 mm ID 18 FR	N - QDE	N - QCE	
5.0 mm ID 20 FR	N - QDF	N - QCF	
5.5 mm ID 22 FR	N - QDG	N - QCG	
6.0 mm ID 24 FR	N - QDH	N - QCH	
6.5 mm ID 26 FR		N - QCI	
7.0 mm ID 28 FR	N - QDJ	N - QCJ	N - Q5Z
7.5 mm ID 30 FR	N - QDK	N - QCK	N - Q5I
8.0 mm ID 32 FR		N - QCL	N - Q52
8.5 mm ID 34 FR		N - QDM	
9.0 mm ID 36 FR		N - QCN	
9.5 mm ID 38 FR		N - QCP	
10.0 mm ID 40 FR		N - QCQ	