


### FOLEY CATHETER



Great comfort for urinary drainage and irrigation

Wide range of catheters for therapeutic drainage, suction, infusion

- Latex coated in silicone, for smooth insertion and comfort for user
- 100% Silicone Foley catheter available, all sizes and gauges
- Reinforced junction for greater security
- 6,8,10FR with rigid catheter guide
- Our 3-way catheters are excellent for bladder irrigation
- Fits luer slip and luer lock connectors, all collection and catheter tip syringe systems
- Standard length 40cm, 2 side openings at balloon tip
- Color coded inflating valve
- Symmetric balloon 3-50ml volumes available
- Blister sterile packaging with inner liner
- 5 year expiry date

LATEX REFERENCE	SILICONE REFERENCE	COLOR	GAUGE	DESCRIPTION
N-QOL1	N-Q7H1	BLACK	6FR	2 way, balloon 3-5ml, 30cm
N-QOQ1		LIGHT GRAY	8FR	2 way, balloon 5-10 ml, 30cm
N-QOR3		WHITE	10FR	2 way, balloon 3-5ml, 30cm
N-QOS3		LIGHT GREEN	12FR	2 way, balloon 5-15ml, 40cm
N-QOT3	N-Q7T3	LIGHT GREEN	14FR	2 way, balloon 5-15ml, 40cm
N-QOU3	N-Q7U3	ORANGE	16FR	2 way, balloon 5-15ml, 40cm
N-QOV3		RED	18FR	2 way, balloon 5-15ml, 40cm
N-QOV3		YELLOW	20FR	2 way, balloon 5-15ml, 40cm
N-QOW3		PURPLE	22FR	2 way, balloon 5-15ml, 40cm
N-QOX3	N-Q8V3	BLUE	24FR	2 way, balloon 5-15ml, 40cm
N-QWU4		RED	18FR	3 way, balloon 5-15ml, 40cm
N-QWU4		RED	18FR	3way, balloon 30ml, 40cm
N-QWV4		RED	18FR	3way, balloon 30ml, 40cm
N-QWV4	N-Q8V4	YELLOW	20FR	3 way, balloon 30ml, 40cm
N-QWW4	N-Q8W4	PURPLE	22FR	3 way, balloon 30ml, 40cm
N-QWX4	N-Q8X4	BLUE	24FR	3 way, balloon 30ml, 40cm