

LARYNGEAL MASK



Short anesthesia alternative

- 100% silicone.
- Resistant, soft, atraumatic ring.
- Balloon capacity up to 60 ml.
- Antireflux valve.
- Inflating tube with marking and pilot balloon.
- Universal connector that fits to any ventilation system.
- Firm connection fitting, no leaks.
- Sterile, 5 year expiry date.

REFERENCE	SIZE	DESCRIPTION
N-RL1	Neonatal	Size 1.0 , 6 ml balloon capacity
N-RLY	Neonatal	Size 1.5 , 10 ml balloon capacity
N-RL2	Pediatric	Size 2.0 , 15 ml balloon capacity
N-RLZ	Pediatric	Size 2.5 , 20 ml balloon capacity
N-RL3	Adult	Size 3.0 , 30 ml balloon capacity
N-RL4	Adult	Size 4.0 , 45 ml balloon capacity
N-RL5	Adult	Size 5.0 , 60 ml balloon capacity