

FEEDING TUBE



- Medical grade PVC, soft and transparent
- Luer slip connector (12FR-22FR conical connector)
- 2 side openings and open tip
- 5 year expiry date

REFERENCE	COLOR	DESCRIPTION
N-QFH	DARK BLUE ●	Feeding tube, 6FR, 40cm
N-QFL	LIGHT BLUE ●	Feeding tube, 8FR, 40cm
N-QFQ	CLEAR ○	Feeding tube, 10FR, 40cm

Allows gastric decompression, sampling, wash, drainage and feeding

- Universal conical connector
- Flexible tube with open tip and 4 side openings
- Markings on tube for accurate insertion
- X- ray line available (REF: N-QI*)
- 5 year expiry date

STOMACH LEVIN TUBE



REFERENCE	COLOR	DESCRIPTION
N-QACH	LIGHT GREEN ●	Stomach Levin tube 6FR X 90cm
N-QHL	BLUE ●	Stomach Levin tube 8FR X 125cm
N-QHQ	BLACK ●	Stomach Levin tube 10 FR X 125cm
N-QHR	WHITE ●	Stomach Levin tube 12 FR X 125cm
N-QHS	GREEN ●	Stomach Levin tube 14 FR X 125cm
N-QHT	ORANGE ●	Stomach Levin tube 16 FR X 125cm
N-QHU	RED ●	Stomach Levin tube 18 FR X 125cm
N-QHV	YELLOW ●	Stomach Levin tube 20 FR X 125cm