


STOPCOCK

Excellent flow control for medicine and biologicals administration



- Available in white, red or blue color
- Rotational knob with flow arrows (360°)
- Connects to all luer slip and luer lock devices
- With extension sets, any length available
- Produced with medical grade PVC, transparent
- High resistance to mechanical or manual flow pressure
- Sterile individual blister packaging
- 5 year expiry date

REFERENCE	DESCRIPTION
N-QTAE	3-way stopcock with extension, 15cm
N-QTAF	3-way stopcock with extension, 50cm
N-QTAG	3-way stopcock with extension, 100cm
N-QUAJA	3-way stopcock, color red
N-QUAJB	3-way stopcock, color blue