


### VAGINAL SPECULUM



- **Siliconized, biocompatible, free of irregular edges**
- **Adjustable press-up trigger**
- **Available in sizes: S, M, L**

#### PRODUCT DESCRIPTION

Size S , 104 x 22 x 8.5mm

Size M , 107 x 31 x 9mm

Size L , 138 x 32 x 9.5mm

**Nubenco Enterprises, Inc.**

Need to contact us? Call 1.201.967.9000 or

Visit our website [www.nubenco.com](http://www.nubenco.com) • e-mail: [info@nubenco.com](mailto:info@nubenco.com)