

Continuous and precise infusion of liquids and parenteral solutions.

MACRODRIP SET



REFERENCE	CHARACTERISTICS
N - AAACQ	Macrodrrip I.V. set 20 drops/ml, bulb, Y-port, luer slip with needle 21G x 1 1/2".
N - ABC - LF	Macrodrrip I.V. set 20 drops/ml, Y-port, luer lock, 1.8cm length, latex-free
N - AA2C	Macrodrrip I.V. set 20 drops/ml, bulb, Y-port luer slip.
N-AAVC	Macrodrrip I.V. set 20 drops/ml, roller and Y clamps, needless Y-port, rotating luer lock, 1.8cm length, DHEP-free
N-AATC-LF	Macrodrrip I.V. set 20 drops/ml, solution filter, Y-port, rotating luer lock (combo lock), 1.8m length, Latex-free
N - AAFC - LF	Macrodrrip I.V. set for pumps 20 drops/ml, LL/LS combo connector. Latex-free
N - AAFE - LFP	Macrodrrip radiopaque I.V. set for pumps 20 drops/ml, LL/LS, combo connector.
N-AB3E-LF	Macrodrrip I.V. set for pumps 20 drops/ml, needless Y-port, rotating luer lock (combo lock), 1.8m length, Latex-free
N-AAVC-LF	Macrodrrip I.V. set 20 drops/ml, needless Y-port, rotating luer lock (combo lock), 1.8m length, Latex-free

Safety and precision in the administration of intravenous fluids

MICRODRIP SET



- Strong and sharp spike, 3 channels.
- 60 drops/ ml, 20 drops/ ml.
- Standard tube length 180 cm.
- Available: Y-port, air inlet, solution filter
- Quantitative roller clamp.
- Sets with luer slip, luer lock or combo (LL/LS) connector.
- Latex-free.
- Easy to open packaging.
- 5 Years expiry date.

REFERENCE	CHARACTERISTICS
N - AA4C	Microdrrip I.V. set, 60 drops/ml, Y-port, bulb, luer slip. 1.8m length
N - AALC	Microdrrip I.V. set, 60 drops/ml, Y-port, rotating luer lock (combo lock), 1.8m length,
N-AB5C-LF	Microdrrip I.V. set 60 drops/ml, Y-port, rotating luer lock (combo lock), 1.8m length, Latex-free