

Nubenco®

Administration of solutions and blood components

Standard syringe 0.5 - 1ml



- 3 pcs syringes: barrel, plunger and rubber tip
- Luer slip connector and detachable needle
- For 0.5ml and 1ml, fixed needle available for high precision of insulin, immunization, BCG and tuberculin dosage
- Color-coded needles
- Sharp 3-bevel, siliconized needle for less traumatic puncture
- Long rubber tip to reduce dead space and residue of medicine
- Siliconized for easy, smooth plunger slide and strong vacuum
- Positive plunger stop
- Permanent engraved fine scale markings on barrel
- Double scale available
- Short scale for BCG to avoid overdose
- Plunger flat top with grids for controlled injection
- Latex or latex-free stopper
- PE envelope or peel-off blister packaging
- 5 year expiry date
- Also available A.D. and safety syringes

REFERENCE	GAUGE	LENGTH	COLOR	APPLICATION
N-S2F5	30G	1/2"	Yellow	For insulin, cosmetic use (Botox)
N-S2FAI	29G	1/2"	Red	
N-S2F4	28G	1/2"	Blue	Allergy test, mesotherapy intradermic use, subcutaneous
N-S2SD5-B-N	27G	3/8"	Light gray	
N-S2SE1-B-N	26G	3/8"	Brown	
N-S2F1	25G	5/8"	Orange	Vaccine Intradermic and subcutaneous use
N-S9ZAZ	25G	5/8"	Beige	
N-S2SV-B-N	23G	1"	Dark Blue	Tuberculin/Immunization
N-S2ST-B-N	22G	1 1/2"	Gray	
N-S2DA1	27G	1/2"		Immunization Insulin
N-S2S1-B-N	25G	5/8"		
N-S9D3	27G(fixed)	1/2"		
N-S9D1	25G(fixed)	5/8"		
N-S2T1-B-NBP	25G	5/8"		
N-S9F5	30G(fixed)	1/2"		
N-S91E3	31G(fixed)	1/2"		

For feeding, irrigation internal wash

- 60ml syringe
- Universal tip fit different gauges tubes or catheters
- Includes separate luer slip adapter
- Siliconized for easy and smooth plunger slide
- Positive plunger stop
- Grid plunger top for controlled irrigation
- 5 year expiry date

Catheter tip syringe 60ml



REFERENCE	ADDITIONAL CHARACTERISTICS
N-SDR-LF	Syringe 60ml, catheter tip, luer slip adapter

Nubenco Enterprises, Inc.

Need to contact us? Call 1.201.967.9000 or

Visit our website www.nubenco.com • e-mail: info@nubenco.com