

### Plaster Of Paris Bandage



SIZE	REFERENCE
3" X 5 yd	N-G3E
4" X 5 yd	N-G3F
5" X 5 yd	N-GEG
6" X 5 yd	N-G3H
8" X 5 yd	N-G316

- Drying time 24 Hours, for greater patient comfort
- Provides a strong and safe immobilization
- High resistance, avoiding damage and increasing patient comfort
- Low plaster losses, for higher resistance and better consistency
- High creaminess, providing a superior finish
- Radiolucent for X-ray imaging and diagnostics
- Air-tight packaging, preserving the product from humidity
- Fiberglass colors available: blue, green, white

### NUCAST® Fiberglass bandage

- POP setting time: 5 - 8 minutes
- Fiberglass setting time: 3 - 5 minutes



SIZE	REFERENCE
2" x 4 yd	N-GSVB
3" x 4 yd	N-GSWB
4" x 4 yd	N-GSXB
5" x 4 yd	N-GSYB
6" x 4 yd	N-GSZB