

BLOOD COLLECTION TUBES



- **Product parts:** plastic or glass tube, rubber stopper, color-coded cap
- **Gamma radiation sterilization**
- **Colors available:** Yellow, red, lavender, blue, black, green, grey

PRODUCT DESCRIPTION

- Gel tube, color: yellow, Additive: gel and clot activator, material: glass/plastic, size: 2ml, 3ml, 5ml, 7ml, Specs: 13x75mm, 13x100mm
- Plain tube, color: red, Additive: no additive/clot activator, material: glass/plastic, size: 2ml-10ml, Specs: 13x75mm, 13x100mm, 16x100mm
- EDTA tube, color: lavender, Additive: EDTAK2/EDTAK3, material: glass/plastic, size: 2ml, 3ml, 4ml, 4.5ml, 5ml, Specs: 13x75mm, 13x100mm
- Coagulation tube, color: blue, Additive: sodium citrate 1:9, material: glass/plastic, size: 1.8ml, 2.7ml, 3.6ml, 4.5ml, Specs: 13x75mm, 13x100mm
- Sodium citrate tube, color: black, Additive: sodium citrate 1:4, material: glass/plastic, size: 1.6ml, 2.4ml, 3.2ml, Specs: 13x75mm, 13x100mm
- Heparin tube, color: green, Additive: sodium heparin, material: glass/plastic, size: 2ml, 3ml, 5ml, Specs: 13x75mm, 13x100mm
- Glucose tube, color: gray, Additive: sodium fluoride/potassium oxalate, material: glass/plastic, size: 2ml, 3ml, 4ml, Specs: 13x75mm, 13x100mm

Nubenco Enterprises, Inc.

Need to contact us? Call 1.201.967.9000 or

Visit our website www.nubenco.com • e-mail: info@nubenco.com