

Maximum precision for dosification and volumetric administration of liquids and medication.

BURETTE



- Single bevel spike for faster flow
- Strong hanger for greater stability
 - Silicone floater available
 - Blue graduated volumetric scale (Available: 100ml, 150ml)
 - Writable white strip on the chamber
- Latex-free option
- Made with recyclable materials
- Detailed instruction on packaging for Complete and incomplete versions

REFERENCE	CHARACTERISTICS
N - AA9F-LF	Burette, macrodrip 20 drops/ml, latex-free
N - AQF	Burette, macrodrip 20 drops/ml, floater
N - AADD	Burette complete set, microdrip 60 drops/ml, LL/LS combined connector
N - AHD	Burette complete set, microdrip 60 drops/ml, Luer slip connector
N-A8B	Burette complete set, microdrip 60 drops/ml, floater, Y-port, Luer lock connector
N-AADD-LF	Burette complete set, microdrip 60 drops/ml, LL/LS combined connector - Latex free