

HYPODERMIC NEEDLE

Easy puncture, minimal pain



REFERENCE	COLOR	DESCRIPTION
N-N1B	White	Needle 16G x ½"
N-N1A	White	Needle 16G x 1"
N-N19	White	Needle 16G x 1½"
N-N1F	Pink	Needle 18G x 1"
N-N1G	Pink	Needle 18G x 1½"
N-N1I	Beige	Needle 19G x 1"
N-N1J	Beige	Needle 19G x 1½"
N-N1L	Yellow	Needle 20G x 1"
N-N1M	Yellow	Needle 20G x 1½"
N-N1P	Green	Needle 21G x 1"
N-N1Q	Green	Needle 21G x 1½"
N-N1S	Black	Needle 22G x 1"
N-N1T	Black	Needle 22G x 1½"
N-N1V	Blue	Needle 23G x 1"
N-N1W	Blue	Needle 23G x 1½"
N-N1Y	Purple	Needle 24G x 1"
N-N1AE	Orange	Needle 25G x 1"
N-N11	Orange	Needle 25G x 5/8"
N-N12	Brown	Needle 26G x ½"
N-N13	Gray	Needle 27G x ½"
N-N1AI	Salmon	Needle 29G x ½"
N-N15	Ice	Needle 30G x ½"
N-N18	Blue	Needle 23G x 1¼"
N-N1X	Purple	Needle 24G x ¾"

- Fits luer slip and luer lock connector
- Highgrade stainless steel, "painless"
- Sharp bevel and siliconized cannula
- Various gauges 16-30G
- Color coded hub
- Barcode printed on individual blister
- 5 year expiry date

Insulin pen needle

- Stainless steel (AISI 304)
- 3 bevelled tip
- Thin wall
- Siliconized
- Compatible with insulin pens from Novo Nordisk and Aventis
- Compatible with EN ISO 11608-2 standard