

URETHRAL NELATON CATHETER



For drainage, urine sampling,
bladder decompression

- Medical grade PVC, soft and transparent
- Universal conical connector (luer slip available for 6, 8, 10FR)
- All gauges available
- 2 side openings
- 5 year expiry date

REFERENCE	PRODUCT	DESCRIPTION
N-QKF	YELLOW	Nelaton catheter, 5 FR, 40cm
N-QKH	LIGHT GREEN	Nelaton catheter, 6 FR, 40cm
N-QKL	BLUE	Nelaton catheter, 8 FR, 40cm
N-QKQ	BLACK	Nelaton catheter, 10 FR, 40cm
N-QKR	WHITE	Nelaton catheter, 12 FR, 40cm
N-QKS	DARK GREEN	Nelaton catheter, 14 FR, 40cm
N-QKT	ORANGE	Nelaton catheter, 16 FR, 40cm
N-QKU	RED	Nelaton catheter, 18 FR, 40cm
N-QKV	YELLOW	Nelaton catheter, 20 FR, 40cm
N-QKW	DARK BLUE	Nelaton catheter, 22 FR, 40cm

Eases controlled suction of
bodily fluids and secretions

- Thumb-controlled connector, fits any suction system
- Easy gauge identification, color-coded
- Gauges: 5,6,8,10,12,14,16, 18 FR x 50 cm
- Medical grade PVC, transparent
- Open tip with 2 side openings
- Markings for controlled insertion
- 5 year expiry date

SUCTION CATHETER



REFERENCE	PRODUCT	DESCRIPTION
N-QEF	YELLOW	Solution catheter, 5 FR, 50cm
N-QEH	GREEN	Solution catheter, 6 FR, 50cm
N-QEL	BLUE	Solution catheter, 8 FR, 50cm
N-QEQ	BLACK	Solution catheter, 10FR, 50cm
N-QER	WHITE	Solution catheter, 12FR, 50cm
N-QES	DARK GREEN	Solution catheter, 14FR, 50cm
N-QET	ORANGE	Solution catheter, 16FR, 50cm
N-QUE	RED	Solution catheter, 18FR, 50cm