

STERILIZATION PAPER



STERIPAP®

Applicable for vapor and ethylene oxide gas. One side medical grade paper and other side transparent film polyester-polyethylene.
Flat or gusseted rolls and pouches in all sizes

Medical Grade Paper

- Easy to use, clear view of products inside
- Longer sterility after sterilization
- 60g/m² good density paper
- Indicators along edges of paper roll
- Repels liquid
- Antimicrobial barrier
- Seals well
- No ETO residue
- Complies with EN 868

REFERENCE	SIZE (metric system)	ROLL TYPE	PACKAGING
N-PF1A	50mm x 200m	FLAT	6 rolls/box
N-PF1B	75mm x 200m	FLAT	4 rolls/box
N-PF1C	100mm x 200m	FLAT	3 rolls/box
N-PF1D	150mm x 200m	FLAT	2 rolls/box
N-PF1E	200mm x 200m	FLAT	1 rolls/box
N-PF1F	250mm x 200m	FLAT	1 rolls/box
N-PF1G	300mm x 200m	FLAT	1 rolls/box
N-PF1H	350mm x 200m	FLAT	1 rolls/box
N-PF1I	400mm x 200m	FLAT	1 rolls/box

REFERENCE	SIZE (metric system)	ROLL TYPE	PACKAGING
N-PG1C	100mm x 100m	GUSSETED	2 rolls/box
N-PG1D	150mm x 100m	GUSSETED	1 rolls/box
N-PG1E	200mm x 100m	GUSSETED	1 rolls/box
N-PG1F	250mm x 100m	GUSSETED	1 rolls/box
N-PG1G	300mm x 100m	GUSSETED	1 rolls/box
N-PG1H	350mm x 100m	GUSSETED	1 rolls/box
N-PG1I	400mm x 100m	GUSSETED	1 rolls/box

REFERENCE	PRODUCT
N-PSP	Sterilization pouch, 40cm x 50cm
N-PSY	Sterilization pouch, 10cm x 20cm
N-PSZ	Sterilization pouch, 20cm x 35cm