

SYRINGES 1 ml - 60ml



REFERENCE	DESCRIPTION
N - S4H	3 mL, Luer Slip, no needle
N - S1HQ	3 mL, Luer Lock, needle 21G x 1½"
N - S1HT	3 mL, Luer Lock, needle 22G x 1"
N - S1HV	3 mL, Luer Lock, needle 23G x 1½"
N - S4I	5 mL, Luer Slip, no needle
N - S1IW	5 mL, Luer Lock, needle 23G x 1½"
N - S2IT-B	5 mL, Luer Slip, needle 22G x 1 ½"
N - S4J	10 mL, Luer Slip, no needle
N - S1JQ	10 mL, Luer Lock, needle 21G x 1½"
N - S1KQ	20 mL, Luer Lock, needle 21G x 1½"
N - S4K	20 mL, Luer Slip, no needle
N - S3M	30 mL, Luer Lock, no needle
N - S3P	50 mL, Luer Lock, no needle
N - S3R	60 mL, Luer Lock, no needle
N - S4R	60 mL, Luer Slip, no needle
N - S1HQ-LF	3 mL, Luer Lock, needle 21G x 1 ½"
N - S2HW	3 mL, Luer Slip, needle 23G x 1 ½"
N - S1IT	5 mL, Luer Lock, needle 22G x 1 ½"
N - S1IQ-LF	5 mL, Luer Lock, needle 21G x 1 ½"
N - S1IT-LF	5 mL, Luer Lock, needle 22G x 1 ½"
N - S1JT	10 mL, Luer Lock, needle 22G x 1 ½"
N - S3J	10 mL, Luer Lock, no needle
N - S1KQ-LF	20 mL, Luer Lock, needle 21G x 1 ½"
N - S3K	20 mL, Luer Lock, no needle
N - S2ST	1mL, Luer slip, 22G x 1 ½"
N - S2SW	1mL, Luer slip, 23G x 1 ½"
N - S2S1	1mL, Luer slip, 25G x 5/8"
N - S2E1	1mL, Insulin 50IU, Luer slip, 26G x 3/8"
N - S2DS	1mL, Tuberculin, Luer slip, 22G x 1 ½"
N - S2DA1	1mL, Tuberculin, Luer slip, 26G x 5/8"

Engineered accuracy in dosage and safety

- 3pcs: Barrel, plunger, rubber tip
- Color-coded attached needle and latex-free optional
- Luer lock or luer slip (centric or eccentric).
- Siliconized for easy and smooth plunger slide.
- Positive plunger stop.
- Permanent engraved fine scale marking on barrel.
- Plunger flat tip with grids for controlled injection
- Breakable plunger - NO REUSE.
- 5 years expiry date.