

Maximum precision for dosification and volumetric administration of liquids and medication.

TRANSFUSION SET



- 15 and 20 drops/ml available.
- Double drip chamber
 - Top chamber with a 200 ug microfilter.
 - 180 cm flexible tube, little memory.
- Calibrated roller clamp for precise drip.
- Latex-free.
- Luer slip connector.
- Easy to open packaging.
- 5 years expiry date.

REFERENCE	CHARACTERISTICS
N - AACB - LF	Transfusion set, 20 drops/ml , luer slip, latex - free.
N - AAIGG	Transfusion set, 20 drops/ml , Y - Port, luer slip, needle 18G x 1 1/2".
N - AJB	Transfusion set, 15 drops/ml , bulb, luer slip.