

### NEBULIZER MASK



Device with optimum mist action for respiratory procedures

- Fitted mask, made with soft PVC
  - Large lateral openings for CO<sub>2</sub> release
  - Adjustable elastic band, latex-free
- Nebulizer cup, 8 ml capacity
- 2m tube for better patient mobility

REFERENCE	ADDITIONAL CHARACTERISTICS
N-Q46	Nebulizer mask, size XL
N-Q47	Nebulizer mask, size L
N-QAK	Nebulizer mask, size M
N-Q49	Nebulizer mask, size S

### VENTURI MASK



- With 2 or 6 color-coded diluters to vary the O<sub>2</sub> concentration.
- With lateral openings on mask for exhaled gas exit
- Designed for maximum comfort and utility
- Latex-free
- Medical grade PVC

REF	DESCRIPTION
N-RV7-6	Adult, 6 diluters
N-RV9-6	Pediatric, 6 diluters
N-RV7	Adult, 2 diluters
N-RV9	Pediatric, 2 diluters

24%	28%	31%	35%	40%	50%
4LPM	4LPM	6LPM	8LPM	8LPM	10LPM