

Nubenco®

Improving care. Preserving life.®

Medical surgical line



- Over 30 years serving the global marketplace
- Produces a wide-range of premium hospital disposables
- All products are developed in the U.S.A. and bear the CE mark
- Many products approved by the US FDA
- Work with Suppliers, Distributors, NGOs, major Multinational Organizations, and Government Agencies worldwide
- Services include product registration, marketing support, training and development, digital media support, and daily one-on-one customer support
- Visit www.nubenco.com to learn about us and how Nubenco can help you become a leader in your market.

Nubenco Enterprises, Inc.
1 Kalisa Way, Paramus, N.J. 07652 U.S.A.

Developed in the U.S.A.

Syringe 200 ml For Contrast Media

Nubenco can supply a wide range of contrast media syringes, to fit all common injectors. We also supply the connectors, both low-pressure connector tubes, and anesthesia extension tubes which are sometimes used as a cheaper substitute for the official connectors. See enclosed. We can even supply the separate filling tube! Lastly, we can supply complete kits with syringe, connector and filling tube.

In some countries, one syringe is used all day for all the patients, changing only the connector and filling tube. In some countries, they even wash, sterilize and reuse the syringe for many days. In some countries, and/or for some low volume contrast media, they do not even use injectors, but administer the media manually through anesthesia extension tubes. So, the business with connectors/filling tubes might be bigger than for the special syringes.

In some parts of the world, the big multinationals control this business by offering "package deals" of some combination of "Free" CT scanner, contrast media, injector and syringe kits. However, even under those conditions, there may be opportunities to supply just some connector and filling tubes, once the customer complies with the minimum required by his package deal. Please, therefore, study your market and advise which injectors are used, whether the business for the consumables is limited due to package deals, and the resultant potential for our syringe kits.

Syringe 200 ml For Contrast Media



- High quality Polypropylene material
- For administration of Contrast Media.
- Kit for CT injector systems:
 - Syringe 200 ml
 - Low pressure connector tube
 - Filling tube
 - Detailed insert
 - Individually boxed
- Fits perfectly to CT injectors such as Medrad 200/300/Vistron Series and Stellant Models
- Kits also available for other systems
- Connector tube with female/male adapter; 150 cm Length available separately (Ref. N-Q2)
- Double packaging for better protection and storage.
- Stands pressures up to 355 PSI/2446 kPa
- 5 year expiry date

REFERENCE	DESCRIPTION
N - SCY	200 ml Syringe, quick fill tube, low pressure tube (fits Medrad 200/300 series)
N - SEY	200 ml Syringe, quick fill tube, low pressure tube (fits Stellant Type Injector)
N - SFY	200 ml Syringe, Frontal injector, quick fill tube, low pressure tube (fits CT9000 Injector)
N - SGY	2-200 ml Syringe, quick fill tube, low pressure tube (fits Stellant Type Injector)

Extension Tubes

NUFLOW flow regulator with extension set



REFERENCE	PRODUCT
N-Q9A-LF	Flow regulator with extension set

T-port extension




REFERENCE	PRODUCT
N-Q1	T-port extension set

Anesthesia line



REFERENCE	PRODUCT
N-QAD	Anesthesia line, luer lock/luer slip, adult
N-QAE	Anesthesia line, luer lock/luer slip, pediatric

Contrast media connectors



REFERENCE	PRODUCT
N-Q6	Filling tube (quick and manual)
N-Q2	Low pressure tube, spiral

Nubenco®

Improving care. Preserving life.®

Thank you!
Gracias!
Merci!
спасибо!
شكرا

Nubenco Group Media – Plug in, Connect, and Join us!

- Visit us sign up for our newsletter! ● www.nubenco.com
- Join us on Facebook! ● <https://www.facebook.com/nubenco>
- Follow us on Twitter! ● <https://twitter.com/NubencoGroup>
- Connect with us on LinkedIn ● <http://www.linkedin.com/company/nubenco-enterprises-inc->
- Subscribe to us on YouTube ● <http://www.youtube.com/user/nubencogroup>

Nubenco Enterprises, Inc.

Need to contact us? Call 1.201.967.9000 or

Visit our website www.nubenco.com ● e-mail: info@nubenco.com